

せんねん灸

SENNENKYU NO KISEKI

Reduced Smoke Moxibustion Products

Instructions

Please be sure to read these instructions before use.

This document includes prohibitions, precautions for use, storage and handling precautions, and usage instructions.

Please read this document carefully and use the product correctly.

Please keep this document in safe place so that it can be easily referred to if necessary.

! Important Precautions

- Moxibustion may result in burns (blister).
- It may also leave marks.

Prohibitions

1. Do not use if you are:

- Persons who are unable to remove the product from the skin themselves
- Children and infants


2. Do not use when you are:

- Suffering from a fever Sweating
- Taking bath 30 minutes before and after moxibustion

3. Do not use on the following locations:

- Face • Mucous membranes
- Areas of eczema, skin irritation, or open wounds
- Newly injured parts of body, including bruises and sprains
- Skin that is wet for reasons such as sweat

4. Do not use in the following manner:

- Do not place the lit end down. Smoke and heat may cause burns (blister) 
- Do not place products on many different locations at the same time.
- Do not apply the product to the same spot repeatedly.
- Do not put foreign objects inside the base of the product.

5. Do not reuse the product.

The product is made of carbonized moxa, so it produces approximately 20% or less of the smoke produced by our other products (third-party tested). For this reason, we recommend it to those who are concerned about smoke.

SENNENKYU no Kiseki is made of carbonized moxa.

Precautions for Use

The following persons should consult a physician or pharmacist prior to use:

- Persons who have had an allergic reaction to a medicine or cosmetic product, including rashes, redness, itchiness and irritation
- Persons whose skin becomes easily irritated due to a skin disorder, sensitivity or imbalance.
- Persons who have a respiratory disease
- Pregnant women
- Persons who have circulatory disorders or who have limited sensitivity to heat, such as diabetics

Cautions During Use

- If it feels hot, remove the product immediately. If it is not removed swiftly, it may result in burns (blister) or leave marks.
- People's perception of heat varies, and it is influenced by product usage conditions, which includes factors such as the environment, the person's tolerance and physical condition, and the location the product was applied.
- Changes in usage conditions, including factors such as the environment, the person's tolerance and physical condition, and the location the product was applied, may cause burns (blister) or leave marks.
- When the moxa is burning, the lit end, the base and the smoke are all hot and may result in burns (blister) or leave marks.
- They may result in burns (blister) or leave marks in areas where the skin is sensitive, including the abdomen, thighs and calves.
- Be careful with the ash that may fall when the product is burning or after it has burned. Wait until the base has cooled before putting the ash in an ashtray and removing the product from the skin.
- Your hands may become dirty because the product is made of carbonized moxa.

Storage and Handling Precautions

- Avoid storing in high humidity and do not allow the product to become wet.
- Store with the provided desiccant.
- Store out of the reach of children and infants.
- These products are not edible and should only be used for the purpose of moxibustion.
- Take extra care when lighting the product to prevent fires.
- After use, make sure that the fire is completely extinguished before discarding the product as general waste.

Moxibustion with 100% Natural Mugwort



The moxibustion product is made of 100% natural mugwort.


Moxa has been used since ancient times, and knowledge about using the gifts of nature has been passed down through generations.

SENNENKYU was developed to introduce a moxibustion product that is easier and more pleasant to use.

SENNENKYU no Kiseki uses carbonized moxa, so it produces less smoke.

Usage Instructions

- Observe your skin condition carefully to ensure safety
- When you use the product for the first time, start with a single point using just one moxibustion product.
- It is recommended to apply the product to one, two or three points on your body in a single moxibustion session.
- People's perception of heat varies.
- The product should feel warm on the skin.
- **Please handle with care as carbonized moxa is easy to break.**



How to Choose the Right Moxibustion Product

Start using SENNENKYU no KISEKI Soft at thermal level 2. SENNENKYU uses five thermal levels.

Please choose the level that best suits your needs or occasion.

SOFT ← warming chart → HARD



SENNENKYU NO KISEKI Soft



SENNENKYU NO KISEKI Regular




SENNENKYU NO KISEKI Hard



1

Peel off the protective backing of the base.



2

Attach the base to the fingertip and light the carbonized moxa.

(It may take some time for it to start burning.)

Note: For safety reasons, we recommend that you attach SENNENKYU to a fingertip and bring it toward the flame.




3

Once lit, apply it to the skin.



4

It will continue to produce heat even after the smoke has disappeared.



5

• **If it feels hot:**

If the product starts to feel hot, remove it immediately.

If it is not removed swiftly, it may result in burns (blister) or leave marks.

• **If it doesn't feel hot:**

Touch the base to make sure it has cooled down before removing it from the skin.



Manufacturer and distributor:

SENNENKYU CO., LTD.

Seller:

SENEFA CORPORATION



<https://sennenq.co.jp>

77 UCHIBO-CHO, NAGAHAMA-SHI, SHIGA #526-0244, JAPAN